


Cold, heat, electrical stimulation, and soft-tissue mobilization, all in...

One Tool.





Dynatron Solaris® Plus Series

The Solaris Plus Series offers the most options in any one device by including Ultrasound, 7 Stim waveforms, and the option of adding Tri-Wave Light Therapy and ThermoStim. Solaris Plus is covered by a 2-Year Warranty and backed by the best customer service in the industry.

709*Plus*

5 Channel Combo Unit (D719)

708Plus

3 Channel Combo Unit (D718)

706*Plus*

5 Channel Stim Unit (D716)

705Plus


3 Channel Stim Unit (D715)

Solaris Plus Cart

The Solaris Plus Cart provides ample storage space with 3 large drawers. Constructed with a low center of gravity to prevent tipping, the Solaris Plus Cart is sturdy and aesthetically pleasing.



7030 Park Centre Drive | Salt Lake City, UT | 84121
toll-free [800] 874-6251 | local [801] 568-7000 | fax [801] 568-7711
www.dynatronics.com | email: info@dynatron.com





Power and Flexibility

Solaris Plus delivers 7 Stim waveforms including IFC, Premod, Biphasic, Russian, High Volt, Microcurrent, Direct Current, and Fixed Frequency IFC/Premod.

In addition, Solaris Plus offers 3-Frequency Ultrasound and the new Tri-Wave Light Therapy. Best of all, Solaris Plus has the power to deliver up to 5 channels of Stim, Ultrasound, and Tri-Wave Light Therapy through both Probe and Pads—all at the same time.

Introducing Tri-Wave Light





The FUTURE is NOW!

The first of its kind, the revolutionary Dynatron ThermoStim™ Probe* combines cold or heat with electrical stimulation and soft-tissue mobilization.

Treat a wide variety of conditions, save time,

and enhance patient outcomes.

ThermoStim, the future of therapeutic modalities.

3-Frequency Ultrasound

The only Ultrasound devices in the industry offering 1, 2, and 3 MHz frequencies for the greatest flexibility in depth of treatment.






THERMOSTIM



Target & Target Sweep

By focusing the treatment precisely where it is needed, Dynatronics' unique TARGET feature provides a more effective interferential treatment. Simply glide your finger across the TARGET



touch pad to move the center of interference to the site of your patient's pain.

Impressively Simple

Even with nearly unlimited flexibility, setting up a treatment has never been easier. Some treatments can be set up in as few as two key presses. Treatment parameters are highly visible, easy to read, and changed at the touch of a fingertip.

*Patent Pending